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
Dodatna literatura:

1. Lippi G, Cadamuro J, von Meyer A, Simundic AM; European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Working Group for Preanalytical Phase (WG-PRE). Practical recommendations for managing hemolyzed samples in clinical chemistry testing. Clin Chem Lab Med. 2018;56(5):718-727.
2. Lippi G, Cadamuro J, von Meyer A, Simundic AM; European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Working Group for Preanalytical Phase (WG-PRE). Local quality assurance of serum or plasma (HIL) indices. Clin Biochem 2018;54:112-118
3. Simundic AM, Baird G, Cadamuro J, Costelloe J, Lippi G. Managing hemolyzed samples in clinical laboratories. Crit Rev Clin Lab Sci. 2020;57(1):1-21.
4. Marques-Garcia F. Methods for Hemolysis Interference Study in Laboratory Medicine – A Critical review. EJIFCC. 2020;31(1):85-97.

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5. Cadamuro J, Šimundić AM, Ajzner E, Sandberg S. A pragmatic approach to sample acceptance. *Clin Biochem* 2017; 50: 579-581.
6. CLSI. Hemolysis, icterus and Lipemia/Turbidity Indices as Indicators of Interference in Clinical Laboratory Analysis – Approved Guideline. CLSI document C56-A. Wayne (PA): Clinical and Laboratory Standards Institute; 2012.
7. Thomas L. Haemolysis as influence and interference factor. *EJIFCC*. 2002;13(4):95-98.
8. Luksic AH, Gabaj NN, Miler M, Dukic L, Bakliza A, Simundic AM. Visual assessment of hemolysis affects patient safety. *Clin Chem Lab Med*. 2018; 56(4): 574-581.
9. Gidske G, Aakre KM, Rustad P, Sandberg S, Norling A, et al. Handling of hemolyzed serum samples in clinical chemistry laboratories: the Nordic hemolysis project. *Clin Chem Lab Med*. 2019;57(11):1699-1711.
10. Šimundić AM, Nikolac N, Ivankovic V, Ferenec-Ružić D, Magdić B, Kvaternik M, Topič E. Comparison of visual vs. Automated detection of lipemic, icteric and hemolyzed samples: can we rely on a human eye? *Clin Chem Lab Med*. 2009; 47(11):1361-1365.
11. Von Meyer A, Cadamuro J, Lippi J, Šimundić AM. Call for more transparency in manufacturers declarations on serum indices: On behalf of the working group for Preanalytical Phase (WG-PRE), European Federation Of Clinical Chemistry and Laboratory Medicine (EFLM). *Clin Chim Acta*, 2018; 484:328-332.